Autor / Author

Unternehmer /
entrepreneur

Status / Status

Aktuell / newsworthy

Kategorie / Category

Produktbeschreibung /
product description

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Economically useful CAD/CAM solution for practice and laboratory

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By now, the dental industry offers a variety of different CAD/CAM systems. For their investment decision, users interested in purchasing a system should consider several criteria. Some of the relevant factors to be considered are the costs, the volume of the scan area, the user friendliness of the design software, the possibility to virtually analyze the occlusion, the range of indications as well as the available selection of material.

Especially the investment in a production unit is cost-intensive – regardless of the manufacturer the user decides on. Generally speaking: the higher the quality requirements are, the bigger the capital invested is. For the average dental practice, optimized systems, which often cost 100,000 Euro and more, are unaffordable. Therefore, many companies offer dentists who would like to integrate technological advancement into their practice or practice laboratory at a reasonable financial risk an adequate option: they can purchase CAD systems consisting of a model scanner and the software but without the milling machine. The scanning of the models and the digital design of the restorations can be done immediately in the dental practice. Afterwards, the data is sent to a milling center for manufacturing. Depending on the material in use, the full-contour monolithic restorations or the anatomically reduced frameworks are sent to the dental practice after a short period of time. There, they only need to be fitted onto the model or processed further according to individual requirements. Dentures milled from material that is adjusted to the color of the residual dentition, can be finished in a simple, efficient and, therefore,

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cost-effective way in the practice laboratory. With challenging patients, alternatively, milled frameworks can be individualized and characterized in the conventional way.

This concept of partial outsourcing is interesting from an economical point of view, since each work step that is done in-house in the practice laboratory presents an additional source of income.

absolute Ceramics matchpoint

absolute Ceramics (biodentis, D-Leipzig) offers dentists and dental technicians with its new scan and design system absolute Ceramics matchpoint the opportunity to partly shift the production of prosthetics into their dental practice without having to invest in expensive complete systems. The stripe-light scanner (Fig. 1), which is part of the system, is capable of digitizing cast models and tooth stumps as well as e.g. bite registrations and partial or circular impressions. Every material which is non-reflective or absorbable – except for black or transparent material – can be captured.

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Fig 1: Stripe light scanner of the scan and design system absolute Ceramics matchpoint.

Along with the scan software, a user-friendly design software is available. With itself, inlays, onlays, partial crowns, caps, crown copings and bridge frameworks, full-contour crowns as well as infix

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Fig. 2: Estimate with the individual logo of the practice laboratory.

crowns and bridges can be designed in a simple and quick way on the provided computer with a convenient 3D mouse.

With a price of 12,800 Euro plus VAT for the stripe-light scanner, scan and design software, computer, monitor, keyboard and 3D mouse, absolute Ceramics matchpoint is suitable for dentists who focus on efficient work processes. Attractive process steps, which have previously been performed by third party laboratories, are now realized in the dental practice. The costly and time-consuming milling or grinding is done in an industrial milling center. In this way, important technical work steps like the modelation and individualization take place in the practice or practice laboratory. This approach is indeed reasonable, since the relevant processes are implemented in direct proximity of the patient as well as the dentist, who can judge best the patient's desires and hopes with respect to the restoration.

To ensure a technically correct handling of the CAD system, an intensive training for two participants is included in the purchase price. In case users have any questions at a later stage, they are offered competent support by the expert team at absolute Ceramics.

Usage

For the arrangement of the scan and modeling process with absolute Ceramics matchpoint, all relevant data is entered into an order form – the so-called order manager. For example, this includes the patient name, the tooth which needs treatment, the kind of restoration, the material, the color system and the tooth shade. On the basis of this information, an estimate with the specific logo of the practice laboratory can be saved as a PDF file and printed out as well. This estimate can be integrated into the full settlement after having finished the order (Fig. 2). Furthermore, it is very convenient that the module for estimates facilitates the overview of outstanding orders from the dental laboratory in order to follow up patients.

Scanning

After the impression, bite-taking and model production, absolute Ceramics matchpoint is used for

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scanning and modeling. When developing the system, special efforts have been made to reduce the manual work steps to a minimum to ensure maximum user comfort.

The scanning process is done in three aligned intervals: overview scan of the cast model, bite registration scan and scan of the tooth stumps. In order to scan the model from different viewing angles for the overview scan, the scan plate implements all relevant movements on its own with the fixed object in the different spatial axes inside the scan chamber. The optical measuring head works on the principle of the linear projection procedure with an accuracy of 12 to 15 µm. At this, stripe patterns are projected onto the surface of the object under measurement. The direction and width of the stripes are dependent on the shape of the surface. A laterally positioned camera in the scan chamber documents these stripes and the system generates the surface shapes with an edge definition of 75 µm according to the stripe pattern. Afterwards, all

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single measurements are automatically composed to a complete data set by the software. During the scanning procedure, 500,000 to 750,000 pixels per minute are recorded in order to construct a virtual copy of the model which is as precise as possible.

For the representation of the opposing dentition and the static occlusion, a bite registration is captured afterwards. The scanning process is completed with the digitization of the free prepared single tooth stumps because the interdendum as well as the preparation margin of tooth stumps which stand closely together are often not clearly visible.

Afterwards, the data of the three single measurements is automatically recorded in one data set and the three-dimensional image of the object can be directly rotated and viewed on the monitor. After having confirmed the scan data they are loaded into the design software.

Modeling

From experience, the design software of the system absolute Ceramics matchpoint enables a specifically intuitive handling. This is due to the illustrative representation of the actual situation in the mouth, the marking of each work step by a self explaining symbol related to the tool in use and the guidance through the complete design process. The workflow is, in principle, similar to the conventional processes in the dental laboratory – merely and simply on a digital way.

The entire modeling process for a crown or partial crown takes an average of not more than ten minutes which is, among other things, due to different automatic modes. For example, the system autonomously determines the preparation margin which can be modified individually using the fine-tuning-mode (Fig. 3). Parameters such as the cement line as well as the correction of the milling radius are predefined but can still be readjusted manually if necessary (Fig. 4) and undercuts are automatically blocked out through a digital procedure. Many different tooth forms are available in a database for a virtual wax-up and can be selected in consideration of the opposing bite (Fig. 5 and 6). Occlusion and interproximal contacts can be specified with a digital wax knife and uneven areas compensated

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Fig. 3: After the software has generated the most probable course of the finishing line, it can be modified with the fine-tuning-mode.

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Fig. 4: Predefined parameters e.g. for the cement line as well as the correction of the millig radius can be readjusted manually.

with the function of smoothing (Fig. 7). If necessary, several restorations can be designed in a 3D view at

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Fig 5: In consideration of the opposing bite, the suitable tooth shape can be selected from a database...

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Fig. 7: With the digital wax knife and the smoothing mode, the occlusal as well as interproximal contacts can be defined more precisely and unevenness is compensated.

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Fig. 8: With the infix technology the crown coping or bridge framework is designed anatomically to ideally support the porcelain work.

the same time so that it is possible to systematically shift the interproximal contacts at a later stage. When the restoration is finished, the data only needs to be exported and with the help of the order manager is then sent to the production center via data stream.

Advantages of the centralized production

Apart from the economical advantages of high production volumes, centralized production centers usually offer a high-precision computer-aided manufacture and a high material quality since the material is based on industrially prefabricated blanks. In terms of milling, innovative as well as conventional materials can be processed. Apart from zirconia and silicate ceramics, the range of materials by absolute Ceramics includes titanium which is considered the most bio-compatible of all metals.

Especially for the manufacture of titanium frameworks, all advantages of CAD/CAM technology become important: as blocks produced in an industrially standardized process can be used, more homogenous and bubble-free dental frameworks are obtained than the conventional casting technique permits.

Current Issues

For infix restorations the framework structure as well as the veneering layer is produced with CAD/CAM technology and assembled in a sintering process with the help of a glass solder (Fig. 8 to 15). An in-vitro

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study led by Prof. Dr. Joachim Tinschert at the Department of Prosthodontics and Dental Materials at RWTH Aachen University, Germany, already stated that, regarding all-ceramic infix crowns, the fracture of the framework as well as chipping can be significantly reduced by using this method of veneering. Since July 2010 the infix technology is not only released for the veneering of zirconia frameworks with silicate ceramics but also of titanium works with feldspar ceramics.

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Fig. 11: ... and the porcelain work are milled separately.

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Fig. 13: ... are assembled in a sintering process with the help of a glass solder.

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Fig. 14: Veneering layer of glass ceramic (l.) and framework structure of zirconia (r.).

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Fig. 15: The inner contour of the porcelain work fits exactly onto the outer contour of the framework structure.

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- 1995-2001 studies of dentistry at the Charité University Hospital / Humboldt-University Berlin, Germany
- 2002-2004 assistant in Berlin
- 2004-2007 freelance dental work in Berlin
- 2005-2007 responsible manager at two practices for homeless people of MUT Gesellschaft für Gesundheit Berlin
- 2006 conferral of a doctorate degree
- 2007-2008 speaker and dental surveyor in the field of performance and health management (PMI) in Munich, Germany
- 2008-2009 head of biodentis Schulungszentrum GmbH in Munich
- since 2010 chief executive of biodentis Schulungszentrum GmbH in Munich

Kontakt / Contact

andre.hutsky@web.de

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- 2008-2009 Leiter der biodentis Schulungszentrum GmbH in München
- seit 2010 Geschäftsführer der biodentis Schulungszentrum GmbH in München